



# Anti-MCART1 (full length) polyclonal antibody (CABT-BL2339)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Full length Human MCART1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Mitochondrion inner membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	1X PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	SLC25A51 (solute carrier family 25, member 51) is a protein-coding gene, and is affiliated with the lncRNA class. An important paralog of this gene is LOC494141.
---------------------	---

# GENE INFORMATION

Entrez Gene ID	<a href="#">92014</a>
Protein Refseq	<a href="#">NP_219480</a>
UniProt ID	<a href="#">Q9H1U9</a>