



Anti-MCART1 (aa 49-79) polyclonal antibody (CABT-BL2340)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Immunogen | Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 49-79 of Human MCART1. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Mouse, Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB |
| Cellular Localization | Mitochondrion inner membrane. |
| Format | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at 4°C (up to 6 months). For long term storage store at -20°C |

BACKGROUND

| | |
|---------------------|---|
| Introduction | SLC25A51 (solute carrier family 25, member 51) is a protein-coding gene, and is affiliated with the lncRNA class. An important paralog of this gene is LOC494141. |
|---------------------|---|

GENE INFORMATION

| | |
|----------------|---------------------------|
| Entrez Gene ID | 92014 |
| Protein Refseq | NP_219480 |
| UniProt ID | Q9H1U9 |