



# Anti-MCL1 polyclonal antibody (CABT-BL2356)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide corresponding to amino acids 121-139 of human MCL1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunoaffinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In 10 mM PBS, pH 7.4 (1% BSA, 0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C.

## BACKGROUND

<b>Introduction</b>	This gene encodes an anti-apoptotic protein, which is a member of the Bcl-2 family. Alternative splicing results in multiple transcript variants. The longest gene product (isoform 1) enhances cell survival by inhibiting apoptosis while the alternatively spliced shorter gene products (isoform 2 and isoform 3) promote apoptosis and are death-inducing. [provided by RefSeq, Oct 2010]
---------------------	--

## GENE INFORMATION

Entrez Gene ID [4170](#)

---

Protein Refseq [NP\\_001184249](#)

---

UniProt ID [A0A087WT64](#)

---