



Anti-MCL1 polyclonal antibody (CABT-BL2358)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from internal region of human MCL-1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Cytoplasm, Nucleus
Preservative	None
Storage	2-8°C

BACKGROUND

Introduction	This gene encodes an anti-apoptotic protein, which is a member of the Bcl-2 family. Alternative splicing results in multiple transcript variants. The longest gene product (isoform 1) enhances cell survival by inhibiting apoptosis while the alternatively spliced shorter gene products (isoform 2 and isoform 3) promote apoptosis and are death-inducing. [provided by RefSeq, Oct 2010]
--------------	--

GENE INFORMATION

Entrez Gene ID	4170
----------------	----------------------

Protein Refseq [NP_001184249](#)

UniProt ID [A0A087WT64](#)
