



Anti-MDP1 (full length) polyclonal antibody (CABT-BL2369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length MDP1 protein (Human)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	MDP1 (magnesium-dependent phosphatase 1) is a protein-coding gene. GO annotations related to this gene include protein tyrosine phosphatase activity and metal ion binding. An important paralog of this gene is NEDD8-MDP1. Magnesium-dependent phosphatase which may act as a tyrosine phosphatase (By similarity).
--------------	---

GENE INFORMATION

Entrez Gene ID [145553](#)

Protein Refseq [NP_001186751](#)

UniProt ID [Q86V88](#)
