



# Anti-MEAF6 polyclonal antibody (CABT-BL2371)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC-P
Cellular Localization	Nucleus; nucleolus
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	Chromatin modification-related protein MEAF6 is a protein that in humans is encoded by the MEAF6 gene.
--------------	--

# GENE INFORMATION

Entrez Gene ID [64769](#)

---

Protein Refseq [NP\\_073593](#)

---

UniProt ID [Q9HAF1](#)

---