



Anti-MEAF6 polyclonal antibody (CABT-BL2372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human C1orf149.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Chromatin modification-related protein MEAF6 is a protein that in humans is encoded by the MEAF6 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [64769](#)

Protein Refseq [NP_073593](#)

UniProt ID [Q9HAF1](#)
