



Anti-mec-2 (aa 19-118) polyclonal antibody (CABT-BL2373)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: GSVKVETRLVSNERSSSIQQEGAMLPSSSSKDDDLLSTSS DEVENMATRTLQGLEESTSIISANSDDDSVKKEKQAEKDV EKGNGKEEKANIQNEFGVCG, corresponding to amino acids 19/118 of Caenorhabditis elegans MEC2
Isotype	IgG
Source/Host	Mouse
Species Reactivity	C. elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell Membrane
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [180826](#)

Protein Refseq [NP_741797](#)

UniProt ID [Q27433](#)
