



Anti-mec-5 (aa 103-202) polyclonal antibody (CABT-BL2375)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: GPTGMKGDKGSMGFPGQKGEGGNSGIPGLKGD TGIPGKVG EPGSTGQPGSKGEKGMGLPGTNGLPAGWP GSKGEDGL PGRPGSPGFP GKKGDVGNGG, corresponding to amino acids 103/202 of Caenorhabditis elegans MEC5
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Caenorhabditis elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [181678](#)

Protein Refseq [NP_741941](#)

UniProt ID [G5EFQ9](#)
