



Anti-MFSD1 (internal region) polyclonal antibody (CABT-BL2383)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to internal sequence amino acids of Human MFSD1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC-P
Cellular Localization	Membrane.
Format	Liquid
Buffer	99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	MFSD1 (major facilitator superfamily domain containing 1) is a protein-coding gene. Diseases associated with MFSD1 include salivary gland cancer.
---------------------	---

GENE INFORMATION

Entrez Gene ID [64747](#)

Protein Refseq [NP_001161375](#)

UniProt ID [Q9H3U5](#)
