



Anti-MFSD7B (aa 223-272) polyclonal antibody (CABT-BL2385)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 223-272 (GPKEVSTACA TAVLGNQLGT AVGFLLPPVL VPALGTQNST GLLAHTQNNT) of Mouse FLVCR.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [226844](#)

Protein Refseq [NP_001074728](#)

UniProt ID [B2RXV4](#)
