



# Anti-MGA (aa 2450-2500) polyclonal antibody (CABT-BL2389)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Immunogen corresponds to a region within amino acids 2450-2500 of Human MGA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IP
<b>Cellular Localization</b>	Nuclear
<b>Format</b>	Liquid
<b>Buffer</b>	Tris citrate/phosphate buffer, pH 7-8
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4°C.

## BACKGROUND

<b>Introduction</b>	Megestrol acetate is used mainly as an appetite stimulant in a variety of conditions and as an antineoplastic agent in the treatment of breast, endometrial, and prostate cancers.
---------------------	--

# GENE INFORMATION

Entrez Gene ID [23269](#)

---

Protein Refseq [NP\\_001157745](#)

---

UniProt ID [B9EGR5](#)

---