



Anti-MGA polyclonal antibody (CABT-BL2390)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human MGA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Megestrol acetate is used mainly as an appetite stimulant in a variety of conditions and as an antineoplastic agent in the treatment of breast, endometrial, and prostate cancers.
--------------	--

GENE INFORMATION

Entrez Gene ID [23269](#)

Protein Refseq [NP_001074010](#)

UniProt ID [B9EGR5](#)
