



# Anti-MIA3 (aa 1-120) polyclonal antibody (CABT-BL2401)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein fragment corresponding to a region within amino acids 1 and 120 of MIA3 (XP_496436).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF
<b>Cellular Localization</b>	Endoplasmic reticulum membrane. Localizes at endoplasmic reticulum exit sites. After loading of COL7A1 into transport carriers, it is not incorporated into COPII carriers and remains in the endoplasmic reticulum membrane.
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	40% Glycerol, 1X PBS, pH 7
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [375056](#)

---

Protein Refseq [NP\\_940953](#)

---

UniProt ID [Q5JRA6](#)

---