



## Anti-MIB2 polyclonal antibody (CABT-BL2402)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide, corresponding to a region between residue 925 and 975 of Human Skelethrophin / Mib2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IP
<b>Cellular Localization</b>	Cytoplasmic
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	Tris citrate/phosphate, pH 7 to 8
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4°C.

### BACKGROUND

<b>Introduction</b>	E3 ubiquitin-protein ligase MIB2 is an enzyme that in humans is encoded by the MIB2 gene.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [142678](#)

---

Protein Refseq [NP\\_001164157](#)

---

UniProt ID [Q96AX9](#)

---