



Anti-MICALL2 (C-terminal) polyclonal antibody (CABT-BL2405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A 16 amino acid peptide from near the C terminal residues of human anti-MICALL2 protein, NM_024723.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Cytoplasmic; cytoskeleton
Format	Liquid
Buffer	1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	79778
Protein Refseq	NP_891554

UniProt ID

[Q6UWK3](#)
