



Anti-MIPLE2 (aa 57-156) polyclonal antibody (CABT-BL2413)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Vector coding for a partial recombinant fusion protein, corresponding to internal sequence amino acids 57-156 of Fruit fly (<i>Drosophila melanogaster</i>) CG18321. Target sequence used to make the antibody: AETVASTTNS DTHRTRNNRK PNKPQEHSQQ VGAHKSGSRK LVVAENGGAL
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	44787
Protein Refseq	NP_651996

UniProt ID

[Q9W0S4](#)
