



Anti-MITD1 (aa 71-120) polyclonal antibody (CABT-BL2417)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 71-120 (RAENIKKYLD QEKEDGKYHK QIKIEENATG FSYESLFREY LNETVTEVWI) of Human MITD1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Late endosome membrane; Peripheral membrane protein; Cytoplasmic side
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [129531](#)

Protein Refseq [NP_620153](#)

