



Anti-MKI67IP polyclonal antibody (CABT-BL2423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A mixture of synthetic peptides corresponding to amino acids 144-160 and 276-293 of human MKI67IP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS (0.2% gelatin, 0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene encodes a protein that interacts with the forkhead-associated domain of the Ki-67 antigen. The encoded protein may bind RNA and may play a role in mitosis and cell cycle progression. Multiple pseudogenes exist on chromosomes 5, 10, 12, 15, and 19.
---------------------	---

GENE INFORMATION

Entrez Gene ID [84365](#)

Protein Refseq [NP_115766](#)

UniProt ID [Q9BYG3](#)
