



Anti-MLLT1 (aa 361-410) polyclonal antibody (CABT-BL2436)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: SEDEASFKSE SAQSSPSNSS SSSDSSSDSD FEPSQNHSQG PLRSMVEDLQ, corresponding to internal sequence amino acids 361-410 of Human MLLT1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Protein ENL is a protein that in humans is encoded by the MLLT1 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [4298](#)

Protein Refseq [NP_005925](#)

UniProt ID [Q03111](#)
