



# Anti-MLLT1 (C-terminal) polyclonal antibody (CABT-BL2437)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide from the C-terminal region of human MLLT1 conjugated to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	Protein ENL is a protein that in humans is encoded by the MLLT1 gene.
--------------	---

# GENE INFORMATION

Entrez Gene ID [4298](#)

---

Protein Refseq [NP\\_005925](#)

---

UniProt ID [Q03111](#)

---