



Anti-MLLT1 (aa 325-375) polyclonal antibody (CABT-BL2438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region between amino acids 325-375 of Human MLLT1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	99% Tris buffered saline, 0.1% BSA
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction Protein ENL is a protein that in humans is encoded by the MLLT1 gene.

GENE INFORMATION

Entrez Gene ID [4298](#)

Protein Refseq [NP_005925](#)

UniProt ID [Q03111](#)
