



Anti-MOGS (N-terminal) polyclonal antibody (CABT-BL2492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH: 10-30 amino acids in length at the N terminal of the first 50 amino acids of human GCS1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ICC/IF, ELISA, WB
Cellular Localization	Endoplasmic reticulum membrane; Single-pass type II membrane protein.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Cleaves the distal alpha 1,2-linked glucose residue from the Glc(3)Man(9)GlcNAc(2) oligosaccharide precursor in a highly specific manner.
---------------------	---

GENE INFORMATION

Entrez Gene ID [7841](#)

Protein Refseq [NP_001139630](#)

UniProt ID [Q13724](#)
