



Anti-MORGUE (aa 463-482) polyclonal antibody (CABT-BL2495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: IYEHERERFEQLVRAWTWKY, corresponding to amino acids 463/482 of <i>Drosophila melanogaster</i> Morgue
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Drosophila</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	SCF ubiquitin ligase complex
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [33648](#)

Protein Refseq [NP_608833](#)

UniProt ID [Q9Y165](#)
