



Anti-MORGUE (aa 111-210) polyclonal antibody (CABT-BL2496)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: TENDVDIMGF EDLEDQPAA VPQHPNPQPP MEPQDQQPAQ PVIDLPQPVP EAAAPRRESP ERDFEHERAI NPFLPIKRP RATAPLALPH YMHELADGRA, corresponding to amino acids 111-210 of Fruit fly (<i>Drosophila melanogaster</i>) Morgue
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Drosophila</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	SCF ubiquitin ligase complex
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [33648](#)

Protein Refseq [NP_608833](#)

UniProt ID [Q9Y165](#)
