



Anti-MPDZ (aa 1008-1433) polyclonal antibody (CABT-BL2497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from domain PDZ 6-8 (aa 1008-1433) of human Mupp1 domain protein conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Reconstitution	Reconstitute in 100 µl of sterile water.
Cellular Localization	Plasma membrane
Format	Lyophilised
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Multiple PDZ domain protein is a protein that in humans is encoded by the MPDZ gene.

GENE INFORMATION

Entrez Gene ID	8777
----------------	----------------------

Protein Refseq	NP_003820
----------------	---------------------------

UniProt ID	O75970
------------	------------------------
