



## Anti-MPDZ (aa 775-825) polyclonal antibody (CABT-BL2498)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	Immunogen corresponds to a region within amino acids 775-825 of Human Mupp1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Plasma membrane
Format	Liquid
Size	100 µg
Buffer	Tris citrate/phosphate buffer, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

### BACKGROUND

Introduction	Multiple PDZ domain protein is a protein that in humans is encoded by the MPDZ gene.
--------------	--

# GENE INFORMATION

Entrez Gene ID [8777](#)

---

Protein Refseq [NP\\_003820](#)

---

UniProt ID [O75970](#)

---