



# Anti-MPDZ polyclonal antibody (CABT-BL2500)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------|
| Immunogen          | Recombinant protein corresponding to amino acids of human MPDZ.                                     |
| Isotype            | IgG                                                                                                 |
| Source/Host        | Rabbit                                                                                              |
| Species Reactivity | Human                                                                                               |
| Purification       | Antigen affinity purification                                                                       |
| Conjugate          | Unconjugated                                                                                        |
| Applications       | WB, IHC                                                                                             |
| Format             | Liquid                                                                                              |
| Buffer             | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)                                                   |
| Preservative       | 0.02% Sodium Azide                                                                                  |
| Storage            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|              |                                                                                      |
|--------------|--------------------------------------------------------------------------------------|
| Introduction | Multiple PDZ domain protein is a protein that in humans is encoded by the MPDZ gene. |
|--------------|--------------------------------------------------------------------------------------|

## GENE INFORMATION

|                |                      |
|----------------|----------------------|
| Entrez Gene ID | <a href="#">8777</a> |
|----------------|----------------------|

Protein Refseq [NP\\_003820](#)

---

UniProt ID [O75970](#)

---