



Anti-MPV17L2 (aa 71-120) polyclonal antibody (CABT-BL2504)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 71-120 (MGPFLHYWYL SLDRLFASG LRGFPNVLKK VLVDQLVASP LLGWYFLGL) of human FKSG24
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein (potential)
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [84769](#)

Protein Refseq [NP_116072](#)

UniProt ID [A0A024R7K6](#)
