



Anti-MRFAP1 (aa 1-30) polyclonal antibody (CABT-BL2505)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 1-30 of Human MRFAP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Nucleus. Cytoplasm; perinuclear region. Note: Colocalizes with MORF4L1 to cell nuclei.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	MORF4 family-associated protein 1 is a protein that in humans is encoded by the MRFAP1 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [93621](#)

Protein Refseq [NP_150638](#)

UniProt ID [Q9Y605](#)
