



Anti-Mrgprc (aa 14-28) polyclonal antibody (CABT-BL2506)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 14-28 of rat Mrgprc.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Lyophilized
Buffer	0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ (5 mg BSA, 0.05 mg sodium azide, 0.05 mg Thimerosal)
Preservative	See individual product datasheet
Storage	Store at -20°C on dry atmosphere. After reconstitution with 200 uL of deionized water and concentration will be 500 ug/mL, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	404646
Protein Refseq	NP_001002286.1

