



Anti-MRGPRH polyclonal antibody (CABT-BL2507)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide as a part of rat MRGPRH conjugated to immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	seven transmembrane domain G-protein coupled receptor; most often found in sensory neurons [RGD, Feb 2006]
---------------------	--

GENE INFORMATION

Entrez Gene ID [404641](#)

Protein Refseq [NP_001002281](#)

UniProt ID [Q7TN38](#)
