



Anti-MRPL51 (internal region) polyclonal antibody (CABT-BL2511)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human MRPL51.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	MRPL51 (mitochondrial ribosomal protein L51) is a protein-coding gene. GO annotations related to this gene include protein binding and structural constituent of ribosome.
--------------	--

GENE INFORMATION

Entrez Gene ID	51258
----------------	-----------------------

Protein Refseq [NP_057581](#)

UniProt ID [Q4U2R6](#)
