



Anti-MS4A12 (full length) polyclonal antibody (CABT-BL2513)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length human MS4A12 protein (AAH29793.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The protein encoded by this gene is a cell surface protein found primarily in the apical membrane of colonocytes. Silencing of this gene in colon cancer cells inhibits the proliferation, cell motility, and chemotactic invasion of cells. This gene is part of a cluster of similar genes
---------------------	--

found on chromosome 11. Two transcript variants encoding different isoforms have been found for this gene.

GENE INFORMATION

Entrez Gene ID	54860
----------------	-----------------------

Protein Refseq	NP_001157942
----------------	------------------------------

UniProt ID	Q9NXJ0
------------	------------------------
