



Anti-MTAP2 polyclonal antibody (CABT-BL2515)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from rat or mouse MAP 2a + 2b + 2c conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cytoplasm, cytoskeleton (Probable).
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [17756](#)

Protein Refseq [NP_001035023](#)

UniProt ID [Q3TLQ0](#)
