



Anti-Mtap2 polyclonal antibody (CABT-BL2516)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant mouse MAP 2 (aa 2 - 309).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, ICC, IHC
Preservative	None
Storage	-20°C, avoid repeated freeze/thaw cycles

GENE INFORMATION

Entrez Gene ID	17756
Protein Refseq	NP_001035023
UniProt ID	Q3TLQ0