



Anti-MTMR1 (aa 541-590) polyclonal antibody (CABT-BL2524)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: FLCNCEQQRK KEDVYTKTIS LWSYINSQLD EFSNPFFVNY ENHVLYPVAS, corresponding to C terminal amino acids 541-590 of Human MTMR1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Plasma Membrane; Peripheral membrane protein, cytoplasmic side.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a member of the myotubularin related family of proteins. Members of this family contain the consensus sequence for the active site of protein tyrosine phosphatases.
---------------------	--

Alternatively spliced variants have been described but their biological validity has not been determined. [provided by RefSeq, Jul 2008]

GENE INFORMATION

Entrez Gene ID	8776
----------------	----------------------

Protein Refseq	NP_003819
----------------	---------------------------

UniProt ID	Q13613
------------	------------------------
