



Anti-MUC3B (aa 791-840) polyclonal antibody (CABT-BL2530)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 791-840 (KTTLKEGLQN ASQDANSCQD SQTLCFKPDS IKVNNNSKTE LTPEAICRRA) of Human MUC3B, XP_001125753.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Single-pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [57876](#)