



Anti-MUG1 (native) polyclonal antibody (CABT-BL2531)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native and purified from rat plasma.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction a plasma proteinase inhibitor involved in acute phase inflammatory response [RGD, Feb 2006]

GENE INFORMATION

Entrez Gene ID [497794](#)

Protein Refseq [NP_075591](#)

UniProt ID [Q03626](#)
