



Anti-MUP1 (full length) polyclonal antibody (CABT-BL2532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full-length mouse Mup1 expressed in E.Coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Aliquot and store at -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [17840](#)

Protein Refseq [NP_001156482](#)

UniProt ID [A2CEL1](#)
