



Anti-MXRA7 (aa 35-84) polyclonal antibody (CABT-BL2541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 35 - 84 (PEPARAPPEP APPAEATGAP APSRPCAPEP AASPAGPEEP GEPAGLGELG) of Human MXRA7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [439921](#)

Protein Refseq [NP_001008528](#)

UniProt ID [P84157](#)
