



Anti-MYCBP2 (aa 2425-2475) polyclonal antibody (CABT-BL2545)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region between residues 2425 and 2475 of Human MYCBP2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Size	100 µl
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Probable E3 ubiquitin-protein ligase MYCBP2 also known as myc-binding protein 2 or protein
--------------	--

associates with myc (PAM) is an enzyme that in humans is encoded by the MYCBP2 gene.

GENE INFORMATION

Entrez Gene ID	23077
----------------	-----------------------

Protein Refseq	NP_055872
----------------	---------------------------

UniProt ID	O75592
------------	------------------------
