



# Anti-MYO7B polyclonal antibody (CABT-BL2558)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human MYO7B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

Entrez Gene ID	<a href="#">4648</a>
Protein Refseq	<a href="#">NP_001073996</a>
UniProt ID	<a href="#">Q6PIF6</a>