



## Anti-N monoclonal antibody, clone C034 (CABT-BL2567)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                           |   |
|---------------------------|---|
| <b>Immunogen</b>          | Bovine respiratory syncytial virus (strains 127, SNK and 9007). Human RSV strains A2, Long, Randall and 8/60. |
| <b>Isotype</b>            | IgG1  |
| <b>Source/Host</b>        | Mouse   |
| <b>Species Reactivity</b> | Bovine, Human, Virus  |
| <b>Clone</b>              | C034  |
| <b>Conjugate</b>          | Unconjugated  |
| <b>Applications</b>       | WB, IHC, IP, ELISA, FC  |
| <b>Format</b>             | Liquid  |
| <b>Buffer</b>             | In PBS, pH 7.5 (0.09% sodium azide)   |
| <b>Preservative</b>       | 0.09% Sodium Azide  |
| <b>Storage</b>            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.           |

### GENE INFORMATION

|                       |                             |
|-----------------------|-----------------------------|
| <b>Entrez Gene ID</b> | <a href="#">911726</a>      |
| <b>Protein Refseq</b> | <a href="#">NP_065398.1</a> |

UniProt ID

[Q65474](#)

---