



# Anti-NAT14 (aa 12-41) polyclonal antibody (CABT-BL2574)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	KLH conjugated synthetic peptide designed between amino acids 12-41 from the N-terminal region of Human NAT14
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Membrane; Single pass membrane protein
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C (up to 6 months). For long term storage store at -20°C

## GENE INFORMATION

Entrez Gene ID [57106](#)

Protein Refseq [NP\\_065111](#)

---

UniProt ID [Q8WUY8](#)

---