



Anti-NDD1 (aa 165-264) polyclonal antibody (CABT-BL2588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: LLDSPSTNFFHKTAKTPLRFVTDSNGAQQSTTENPGQ QQNVFSNVDLNNLLKSNKTPSSSCTGAFSRTPLSKIDMN LMFNQPLPTSPSKRFSSLSL, corresponding to amino acids 165/264 of <i>S. cerevisiae</i> Ndd1p.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [854554](#)

Protein Refseq [NP_015017](#)

UniProt ID [Q08887](#)
