



# Anti-NECAB2 (full length) polyclonal antibody (CABT-BL2592)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Full length EFCBP2 protein (Human)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF
<b>Cellular Localization</b>	Cytoplasm.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	N-terminal EF-hand calcium-binding protein 2 is a protein that in humans is encoded by the NECAB2 gene.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [54550](#)

---

Protein Refseq [NP\\_061938](#)

---

UniProt ID [Q7Z6G3](#)

---