



Anti-NEF3 polyclonal antibody (CABT-BL2598)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 800 to the C-terminus of Rat 160 kD Neurofilament Medium.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	neuronal intermediate filament protein; important component of the cytoskeleton [RGD, Feb 2006]
---------------------	---

GENE INFORMATION

Entrez Gene ID [281347](#)

Protein Refseq [NP_001231074](#)
