



Anti-NEURL2 (internal region) polyclonal antibody (CABT-BL2607)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human NEURL2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Neuralized-like protein 2 is a protein that in humans is encoded by the NEURL2 gene.
--------------	--------------------------------------------------------------------------------------

GENE INFORMATION

Entrez Gene ID	140825
----------------	------------------------

Protein Refseq [NP_542787](#)

UniProt ID [Q9BR09](#)
