



Anti-NFXL1 (aa 100-150) polyclonal antibody (CABT-BL2615)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 100-150 of Human NFXL1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Membrane.
Format	Liquid
Size	100 µg
Buffer	99.91% Tris citrate/phosphateNote: pH: 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID [152518](#)

Protein Refseq [NP_694540](#)

UniProt ID [B1Q2B0](#)
